



Manufactured By: HYPHEN BioMed

Cross-Reactivities with Homologous Animal Proteins for the Various ZYMUTEST Kits

Species/Factor	Human	Monkey Cynomolgus	Monkey Rhésus	Pig	Goat	Sheep	Horse	Dog	Rabbit	Guinea Pig	Rat	Mouse
Annexin V	++++	ND	ND	++++	++++	++++	+++	++	-	+++	+	++++
PF4	++++	++++	++++	+/-	++++	+	-	-	-	+++	++	+++
CRP	++++	++++	++++	-	-	-	-	-	++++	-	-	-
tPA Antigen	++++	++++	++	++++	++++	++++	+++	+	++++	+/-	+/-	+/-
PAI-1 Antigen	++++	++	++++	-	+/-	+/-	-	-	-	-	-	-
Protéine S, Free	++++	++++	++++	-	-	-	-	-	-	-	-	-
Protéine S, Total	++++	++	++	-	-	-	-	-	-	-	-	-
β ₂ GPI	++++	++++	++++	+++	+/-	+	++	++	-	+++	+++	++
Fibrinogen	++++	++++	++++	++	+	++	++	+	-	-	+	+
Protein C	++++	++++	++++	+/-	++++	++	-	-	-	+/-	-	-
Fibronectin	++++	+++	++++	++++	++++	++++	++++	++++	-	++++	++++	++++
vWF	++++	++++	++++	++++	++++	++++	++++	++++	-	++	+++	++
Protein Z	++++	++++	++++	+/-	+++	++	+	+	+/-	++	+	++
Factor IX	++++	++++	++++	++++	++++	+++	+++	+++	+/-	++++	++++	++++
Factor X	++++	++++	++++	++	+++	++	+	+	-	+++	++	++
Factor XIII-A	++++	++	++	+++	+++	+	+	+/-	++	-	-	-

Estimation of respective cross-reactivities:

- ++++ > 60 % cross-reactivity
- +++ 40 to 60 % cross-reactivity
- ++ 20 to 40 % cross-reactivity
- + 10 to 20 % cross-reactivity
- +/- 2 to 10 % cross-reactivity
- < 2 % cross-reactivity

Form AH90
08-2008



7768 Service Center Drive • West Chester OH 45069

Phone: 513.770.1991

Toll Free: 866.783.3797

Fax: 513.573.9241

Email: info@aniara.com

www.aniara.com